

## RION Syringe Sampler KZ30W1/30W2

- Choice of three operation modes: purge mode, measurement mode, and combination mode for preset sequence of purge and measurement operations.
- Memory for up to 50 sets of parameter settings for purge mode or measurement mode.
- Selectable blank volume lets sample fluid flow rate stabilize after start of plunger action.
- Syringe set with integrated valve for easy replacement.
- Automatic syringe capacity detection.
- The syringe sampler KZ-30W2 is in conformity with the European Union's CE Compliance requirements.



### Specification

Counters and sensors	KL-11, KL-11A, KL-20A, KL-22, KL-26, KL-27(with some limits regarding usable functions) KS-40/40A/40AF/41, KE-40/40AKE-40 (with some limits regarding usable functions)
Syringe size	25 mL PP syringe set (standard) 25 mL Pylex glass syringe set (option) 10 mL Pylex glass syringe set (option)
Repeat precision	60mm+/-2%
Resolution	0.5% or less
Purge/Measurement volume	Setting range: 10 to 100mL/min (with 25 mL syringe set) 5 to 50mL/min (with 10 mL syringe set) Resolution: 1mL/min
Number of memory addresses:	50 (with battery backup)
Power Supply	AC100-240 V, 50/60 Hz, Approx. 100 VA
Dimensions and Weight	147 (W) X 290 (D) X 350(H) mm Approx. 6 kg

**Your Local Distributor:**



**Please Visit Us At:**

[www.jalttechnologies.com](http://www.jalttechnologies.com)

**P.O. Box 597**

**Manassas, VA 20108**

**703-335-9910**

JALT Technologies is a registered trademark of JALT Technologies, LLC. Copyright 2005. All rights reserved. Rion is a registered trademark of Rion Co. LTD..